## Protocol for annealing oligos

Written by Seok-Yong Choi on 11-28-2008

- 1. Mix 200 ng of each oligo
- 2. Incubate the following reaction for 3 min at 100 °C.

Primer 200 ng Buffer H 5 µl Water up to 50 µl

- 3. Incubate for 5 min at 65-70 °C.
- 4. Cool down the reaction at RT.
- \* In ligation reaction, vector:oligo ratio should be 1:9.

Primer (5  $\mu$ M) 1  $\mu$ l each Buffer H\* 5  $\mu$ l Water 43  $\mu$ l

- $\rightarrow$  3 min at 100 °C.
- $\rightarrow$  5 min at 70 °C.
- → Cool down the reaction at RT
- → Primer concentration: 100 nM

 $20 \text{ ng/}\mu\text{l of vector } (6000 \text{ bp}) = 5 \text{ nM}$ 

\*H buffer: from Roche; restriction enzyme buffer 50 mM Tris HCl 10 mM MgCl<sub>2</sub> 100 mM NaCl 1 mM Dithioerythritol (DTE) pH 7.5 at 37 °C

To convert µg to pmol:

$$\mu g \times \frac{10^6 pg}{1 \mu g} \times \frac{pmol}{330 pg} \times \frac{1}{N} = pmol$$

where N is the number of nucleotides and 330pg/pmol is the average MW of a nucleotide.