

Extraction of genomic DNA from zebrafish embryos for PCR

Written by Seok-Yong Choi on 11-27-2008
Reference: http://zfin.org/zf_info/zfbook/chapt9/9.3.html

1. Add 50 μ l of extraction buffer per embryo.
2. Incubate at 50 °C for 2 - 3 hr or O/N. Tap the tube occasionally during the incubation.
3. Incubate at 95 °C for 10 min to heat-inactivate Proteinase K.
4. Spin at 12,000 rpm for 10 min at 4 °C.
5. Use 2 – 3 μ l to set up PCR reaction in a total volume of 10 μ l.

Optionally, you can use ethanol precipitation to get more pure DNA.

6. Add 100 μ l of EtOH, mix and place on ice for 20-30 min.
7. Centrifuge at full speed at 4 °C for 10 min, remove supernatant and add 200 μ l of 70% EtOH.
8. Spin again for 2 min, remove liquid and air-dry pellet.
9. Resuspend the DNA in 20 μ l of 10 mM Tris (pH 8.0).
10. Proceed with PCR reaction. Usually 1/4 to 1/2 (5-10 μ l) of the DNA isolated from one embryo is sufficient.

Extraction buffer

	Stock	Final	Total 50 ml
Tris HCl (pH 8.0)	1 M	10 mM	0.5 ml
EDTA	125 mM	2 mM	0.8 ml
SDS	10%	0.5%	2.5 ml
Proteinase K*	20 mg/ml	200 µg/ml	0.5 ml

* Add just before use.

	Stock	Final	Total 50 ml
Tris HCl (pH 8.0)	1 M	10 mM	0.5 ml
NaCl		200 mM	
EDTA	125 mM	10 mM	4 ml
SDS	10%	0.5%	2.5 ml
Proteinase K*	20 mg/ml	200 µg/ml	0.5 ml