

# Immunofluorescence Cell Staining

Written by Seok-Yong Choi in 2003

1. Take out the coverslip and drain it onto Kimwipes.
2. Put the coverslip in 6-well plate having 3% paraformaldehyde (PFA) and let stand for 15 min.
3. Wash cells 3X with PBS for 10 min.
4. Replace fixing solution with 0.2% Triton X-100/1% normal goat serum (NGS) in PBS and let stand for 5 min on ice.
5. (Optional) When staining nucleic acid, after 3X 10-min washing with PBS, add dyes, let stand for 5-10 min and wash 2x with PBS, 10 min each
6. Wash cells 3X with 1% NGS in PBS, 10 min each.
7. Transfer coverslip to a piece of parafilm set on a flat surface. Add 100  $\mu$ l of primary antibody diluted in 5% NGS in PBS. Incubate at RT for 1 hr.
8. Remove the antibody solution from a coverslip, and wash 3X with 200  $\mu$ l of 1% NGS in PBS, 10 min each.
9. Add 100  $\mu$ l of secondary antibody diluted in 5% NGS in PBS to coverslips, and incubate at RT for 1 hr.
10. Wash 3X with 200  $\mu$ l of PBS, 10 min each.
11. Mount the coverslip on the slide with Vectashield.
  - 12 mm round coverslip: 3  $\mu$ l of Vectashield or 4  $\mu$ l of Prolong
  - 18 mm rectangle coverslip: 8  $\mu$ l of Vectashield or 10  $\mu$ l of Prolong

\* Approximately 2 ml of 5% NGS in PBS is required for the staining of each coverslip.