

The Inoue Method for Preparation of Competent *E. coli*

Reference: *Molecular Cloning* 1.112

Pick a single bacterial colony (DH5 α) and transfer the colony into 4 ml of LB broth. Incubate the culture overnight at 37°C with vigorous shaking.

Early in the morning, inoculate 125 ml of LB with 0.2, 0.4, 1 or 2 ml of the starter culture. Incubate the culture 24 hr at 18°C with moderate shaking.

The first thing in the following morning is to put a frozen Inoue transformation buffer* (50 ml) on ice in order to thaw it.

When the OD₆₀₀ of the culture reaches 0.55, transfer the Erlenmeyer to an ice-water bath for 10 minutes.

Harvest the cells by centrifugation at 2,500g for 10 minutes at 4°C.

Pour off the medium and store the open centrifuge bottle on a stack of paper towels for 2 minutes. Use a vacuum aspirator to remove any drops of remaining medium adhering to walls of the centrifuge bottle or trapped in its neck.

Resuspend the cells gently in 40 ml of ice-cold Inoue transformation buffer.

Harvest the cells by centrifugation at 2,500g for 10 minutes at 4°C.

Remove the supernatant using the procedure described above.

Resuspend the cells gently in 10 ml of ice-cold Inoue transformation buffer.

Add 750 μ l of DMSO. Mix the bacterial suspension by swirling and then store it in ice for 10 minutes.

Working quickly, dispense aliquots (50 μ l each) of the suspensions into chilled, sterile microfuge tubes (about 200 microfuge tubes are needed). Immediately snap-freeze the competent cells by immersing the tightly closed tubes in a bath of liquid nitrogen. Store the microfuge tubes at -70°C.

* Inoue transformation buffer

PIPES 0.5 M (pH 6.7) (Sigma # P1851, not disodium salt)

PIPES 7.55g in 40 ml of pure H₂O

Adjust the pH to 6.7 with 5 M KOH or HCl.
Add pure H₂O to 50 ml.

Store PIPES in -20 °C.

MnCl ₂ ·4H ₂ O	5.44g
CaCl ₂ ·2H ₂ O	1.10 g
KCl	9.33 g
PIPES (0.5 M, pH 6.7)	10 ml
H ₂ O	to 500 ml

Filter through a prerinsed 0.45-μm Nalgene filter. Store at -20 °C.